

## General Stain-Free Western Blotting Protocol

### Abstract

Western blotting is a widely adopted powerful technique to separate and detect proteins, estimate their quantities, and determine their molecular weights. Proteins of interest are extracted from cells or tissues and processed in a multistep western blotting workflow. This protocol describes Bio-Rad's Stain-Free Western Blotting Workflow, which takes only 5 hr to complete versus 16 hr for traditional western blotting. When following the lengthy and laborious traditional western blotting workflow, researchers assess the adequacy of their western blot data at the final step of the procedure. Conversely, Stain-Free technology enables users to determine transfer efficiency and collect additional data quality information quickly and easily throughout the western blotting process.

### Procedure for Chemiluminescent Western Blotting

#### 1 Protein Sample Preparation for Western Blot Analysis

1. Place the cell culture dish on ice and wash the cells with ice-cold Tris-buffered saline (TBS; 20 mM Tris, pH 7.5, 150 mM NaCl).
2. Aspirate the TBS, then add 1 ml per 100 mm dish ice-cold RIPA buffer (50 mM Tris-HCl [pH 8.0], 150 mM NaCl, 1% NP-40, 0.5% sodium deoxycholate, 0.1% SDS, 1 mM NaF, 1 mM sodium orthovanadate, Roche Protease Inhibitor Tablet).
3. Scrape adherent cells off the dish using a cold plastic cell scraper and gently transfer the cell suspension into a precooled microcentrifuge tube.
4. Maintain constant agitation for 30 min at 4°C.
5. If necessary, sonicate three times for 10–15 sec to complete cell lysis and shear DNA to reduce sample viscosity.
6. Spin at 16,000 x g for 20 min in a 4°C precooled centrifuge.
7. Gently remove the microcentrifuge tube and place it on ice. Transfer the supernatant to a fresh tube, also kept on ice, and discard the pellet.
8. Remove a small volume (10–20 µl) of lysate to perform a protein assay. Determine the protein concentration for each sample.
9. If necessary, aliquot the protein samples for long-term storage at –20°C. Repeated freeze-thaw cycles cause protein degradation and should be avoided.

10. Take 20 µg of each sample and add an equal volume of 2x Laemmli sample buffer.
11. Boil each cell lysate in sample buffer at 95°C for 5 min.
12. Centrifuge at 16,000 x g for 1 min.

#### 2 Electrophoresis with Stain-Free Gels (~30 min)

1. Take an Any kD Criterion TGX Stain-Free Precast Gel, remove the comb, and remove the tape from the bottom of the cassette.
- Note:** Gel percentage selection depends on the size of the protein of interest. A 4–20% gradient gel separates proteins of all sizes and is also a good option.
2. Place the cassette in a Criterion Cell and fill each integrated upper buffer chamber with 60 ml 1x Tris/glycine/SDS running buffer (dilute from 10x TGS, 192 mM glycine, 0.1% SDS, pH 8.3).
  3. Fill each half of the lower buffer tank with 400 ml running buffer to the marked fill line.
  4. Load the protein samples (10–20 µg each) and protein standard.

**Note:** We recommend using Precision Plus Protein Unstained Standards with TGX Stain-Free Gels and Precision Plus Protein All Blue Standards with fluorescent western blots.

5. Place the lid on the tank, aligning the color-coded banana plugs with the corresponding jacks on the lid.
6. Set the voltage to 200–300 V and run the gel for 20–30 min.

**Note:** Midi format gels take approximately 30 min; mini format gels take only 15–20 min.

### 3 Visualize Separation Using the ChemiDoc™ Go or ChemiDoc MP Imaging System (~5 min)

1. After electrophoresis is complete, turn off the power supply and disconnect the electrical leads.
2. Remove the gel cassette from the cell. Pull the two plates of the cassette apart to expose the gel using the lid of the Criterion Cell.
3. Carefully lift the gel from the cassette and place it on the center of the Blot/UV/Stain-Free Sample Tray for the ChemiDoc Go or ChemiDoc MP Imaging System.
4. Capture the Stain-Free gel image using the following settings:
  - Application: Stain-Free Gel
  - Gel activation time: 45 sec (ChemiDoc MP) or blot (ChemiDoc Go)
  - Image Size (W x L): medium; 15.5 x 12.4 cm
  - Exposure: optimal automatic exposure
5. Remove the gel from the UV/Stain-Free Sample Tray and keep it wet in running buffer.

### 4 Protein Transfer with the Trans-Blot™ Turbo System (~10 min)

1. Prepare a membrane from the Trans-Blot Turbo RTA Midi 0.45 µm LF PVDF Transfer Kit or open a Trans-Blot Turbo Midi Nitrocellulose Transfer Pack and place the pad with the membrane on the base of the transfer cassette.
2. Place the gel on top of the membrane, place the top pad on the gel, and roll out bubbles.
3. Place the lid on the cassette base and lock it.
4. Insert the cassette into either instrument bay.
5. Start the transfer by selecting the preset Turbo program and choosing the Criterion Gel size, and then press **RUN**. For most proteins, the 7 min protocol is ideal.
6. When the transfer is complete, disassemble the blotting sandwich and place both the blot and the gel in a container with deionized water.

### 5 Verify Transfer Using the ChemiDoc Go or ChemiDoc MP Imaging System (~5 min)

1. Place the post-transfer gel on the center of the UV/Stain-Free Sample Tray of the ChemiDoc Go or ChemiDoc MP Imaging System.
2. Capture the Stain-Free image of the post-transfer gel with the following settings:
  - Application: Stain-Free Gel
  - Gel activation time: none
  - Image Size (W x L): medium; 15.5 x 12.4 cm
  - Exposure: same as the exposure time for the pretransfer gel image

**Note:** This Stain-Free image of the post-transfer gel is used to verify the transfer efficiency.

3. Remove the gel from the UV/Stain-Free Sample Tray and image the blot with the following settings:
  - Application: Stain-Free Blot
  - Image Size (W x L): medium; 15.5 x 12.4 cm
  - Exposure: optimal automatic exposure

**Note:** Keep the blot wet all the time; do not let it dry. This Stain-Free blot image is used to check the transfer quality.

4. Remove the blotting membrane from the UV/Stain-Free Sample Tray and place it in a container with 1x TBS + 0.1% Tween-20 (TBST).

### 6 Antibody Incubation (~4 hr)

1. Block in EveryBlot Blocking Buffer for 5 min at room temperature with agitation.
2. Dilute the primary antibody in full-strength EveryBlot Blocking Buffer and incubate the membrane for 1 hr at room temperature with agitation.

**Note:** The antibody should be diluted in the blocking buffer according to the manufacturer's recommended ratio. The primary antibody may also be applied to the blot overnight at 4°C, depending on the antibody quality and performance.

3. Rinse the blot five times for 5 min with TBST.
4. Dilute the horseradish peroxidase (HRP)-conjugated secondary antibody in full-strength EveryBlot Blocking Buffer and incubate the membrane for 1 hr at room temperature with agitation.

**Note:** As an alternative to EveryBlot Blocking Buffer, the antibody can be diluted in 1% casein or 5% skim milk in TBS at the manufacturer's recommended ratio. However, blocking times will vary.

5. Rinse the blot five times for 5 min with TBST.

**7 Imaging and Analysis Using the ChemiDoc Go or ChemiDoc MP Imaging System and Image Lab Software (~10 min)**

- After antibody incubation, place the blot back on the UV/Stain-Free Sample Tray of the ChemiDoc Go or ChemiDoc MP Imaging System and capture a Stain-Free image of the blot with the following settings:
    - Application: Stain-Free Blot
    - Image Size (W x L): medium; 15.5 x 12.4 cm
    - Exposure: optimal automatic exposure
- Note:** Keep the blot wet all the time; do not let it dry. This Stain-Free blot image is used for total protein loading control and normalization.
- After imaging, keep the membrane in TBST while preparing the Clarity™ Max Western ECL Substrate mixture.
  - Mix the Clarity Max Substrate Kit components in a 1:1 ratio. Prepare 0.02 ml of solution/cm<sup>2</sup> of membrane. For a midi-sized membrane (8.5 x 13.5 cm), approximately 2 ml of solution is sufficient.
  - Incubate the membrane in the substrate solution for 5 min.
  - Place the blotting membrane back on the UV/Stain-Free Sample Tray of the ChemiDoc Go or ChemiDoc MP Imaging System.
  - Capture the chemiluminescent signals on the blot with the following settings:
    - Application: Chemiluminescent Blot
    - Image Size (W x L): medium; 15.5 x 12.4 cm
    - Exposure: optimal automatic exposure

**Note:** The Image Size for the Stain-Free and chemiluminescence images must be the same for follow-up total protein normalization by Image Lab Software.

**Validate Total Protein Normalization (~5 min)**

- Download the chemiluminescence image of the protein of interest and the Stain-Free image captured in the imaging step and open them in Image Lab Software.
- Link the chemiluminescence and Stain-Free images for normalization analysis by creating a multichannel image.
  - Select **Create Multichannel Image** from the File dropdown menu
  - Drag the chemiluminescence and Stain-Free images into channels 1 and 2

**Note:** The two images may not align perfectly.

- Click **OK**
- Hide the multichannel overlay and chemiluminescence channel by deselecting them from the horizontal toolbar (located at the top of the image)

- Select and use **Detect Lanes and Bands** tools to detect lanes in the Stain-Free image.

**Note:** Although unnecessary for total protein normalization, bands can also be defined if a protein standard is present. Define bands only in the protein standard lane.

- Define the normalization channel.
  - In the normalization section under Lanes and Bands, set the normalization channel to Stain-Free Blot
  - Define a normalization lane

**Note:** This lane should not contain a protein standard.

- Show the chemiluminescence image by selecting it from the horizontal toolbar
  - Select **Sync Normalization Lanes** to populate lanes on the chemiluminescence image and make adjustments if necessary
- Detect the bands of interest in the chemiluminescence image.
  - To view the normalized target protein volumes, click **Analysis Table** on the toolbar. All calculations will be performed automatically by the software, including Normalization Factor and Normalized Volume. The target protein band intensity values are now adjusted for variation in the protein load. This will allow for accurate comparison of target proteins among samples.

**Procedure for Multiplexing Fluorescent Western Blotting**

Follow the steps described in the Procedure for Chemiluminescent Western Blotting section for sample preparation, electrophoresis, protein transfer, transfer verification, and antibody transfer. Variation in protocols between the chemiluminescent and multiplex fluorescent western blotting procedures will start at the secondary antibody incubation step.

**1 Antibody Incubation**

Follow the instructions in steps 1–5 in the Antibody Incubation section on page 2.

- Using full-strength EveryBlot Blocking Buffer, dilute appropriate secondary antibodies (StarBright™ 520 and StarBright 700) corresponding to the species of the primary antibodies and then incubate the membrane for 1 hr at room temperature with agitation.

**Note:** If EveryBlot Blocking Buffer was not used, the antibodies can be diluted in 1% casein or 5% skim milk in TBS according to the manufacturer's recommended ratios. Selection of commercially available secondary antibodies should be based on the species of the primary antibodies and the fluorescence compatibility with the imaging system.

- Rinse the blot five times for 5 min with TBST.

## 2 Imaging and Analysis by Image Lab Software — Total Protein Normalization (~5 min)

- Make the following selections to acquire a three-channel multiplex fluorescence image of the blot.
 

Channel 1:

  - Application: StarBright 520 Blot
  - Image Size (W x L): medium; 15.5 x 12.4 cm
  - Exposure: optimal automatic exposure

Channel 2:

  - Application: StarBright 700 Blot
  - Image Size (W x L): medium; 15.5 x 12.4 cm
  - Exposure: optimal automatic exposure

Channel 3 (optional):

  - Application: Far Red or Near-Infrared (FR/NIR) Blot (DyLight 680 or 800)
  - Image Size (W x L): medium; 15.5 x 12.4 cm
  - Exposure: optimal automatic exposure
- Download the multichannel image captured on the ChemiDoc Go or ChemiDoc MP Imaging System and open it in Image Lab Software. Hide the multichannel overlay and fluorescence channels by deselecting them from the horizontal toolbar (located at the top of the image).
- Select and use **Lane and Bands tools** to detect lanes for the Stain-Free image.

**Note:** Although unnecessary for total protein normalization, bands can also be defined if a protein standard is present. Define bands only in the protein standard lane.

- Define the normalization channel.
  - In the normalization section under Lanes and Bands, set the normalization channel to Stain-Free Blot
  - Define a normalization lane

**Note:** This lane should not contain a protein standard.

- Show the fluorescence channel images by selecting them from the horizontal toolbar
  - Press **Sync Normalization Lanes** to populate lanes on the fluorescence channel images and make adjustments if necessary
- Detect the bands of interest in the fluorescence channel images.
  - To view the normalized target protein volumes, click **Analysis Table** on the toolbar. All calculations will be performed automatically by the software, including Normalization Factor and Normalized Volume. The target protein band intensity values are now adjusted for variation in the amount of protein loaded. This will allow for accurate comparison of target proteins among samples.

Visit [bio-rad.com/StainFreeWestern](http://bio-rad.com/StainFreeWestern) for more information.

BIO-RAD, CHEMIDOC, CLARITY, STARBRIGHT, and TRANS-BLOT are trademarks of Bio-Rad Laboratories, Inc. in certain jurisdictions. All trademarks used herein are the property of their respective owner. © 2025 Bio-Rad Laboratories, Inc.

The ChemiDoc Go Imaging System uses the Linux operating system, which contains software licensed under the following licenses as well as others:

Lesser Gnu Public License v 2.0, 2.1, 3.0

Apache License 2.0

BSD License

Clarity Max Western ECL Substrate and/or its use is covered by the following U.S. patents under license to Bio-Rad Laboratories, Inc., including, but not limited to, U.S. Patent Nos. 7,855,287, 9,040,252, and 8,129,136.



**Bio-Rad  
Laboratories, Inc.**

Life Science  
Group

**Website** [bio-rad.com](http://bio-rad.com) **USA** 1 800 424 6723 **Australia** 61 2 9914 2800 **Austria** 00 800 00 24 67 23 **Belgium** 00 800 00 24 67 23 **Brazil** 4003 0399  
**Canada** 1 800 361 1808 **China** 86 21 6169 8500 **Czech Republic** 00 800 00 24 67 23 **Denmark** 00 800 00 24 67 23 **Finland** 00 800 00 24 67 23  
**France** 00 800 00 24 67 23 **Germany** 00 800 00 24 67 23 **Hong Kong** 852 2789 3300 **Hungary** 00 800 00 24 67 23 **India** 91 124 4029300 **Israel** 0 3 9636050  
**Italy** 00 800 00 24 67 23 **Japan** 81 3 6361 7000 **Korea** 82 080 007 7373 **Luxembourg** 00 800 00 24 67 23 **Mexico** 52 555 488 7670  
**The Netherlands** 00 800 00 24 67 23 **New Zealand** 64 9 415 2280 **Norway** 00 800 00 24 67 23 **Poland** 00 800 00 24 67 23 **Portugal** 00 800 00 24 67 23  
**Russian Federation** 00 800 00 24 67 23 **Singapore** 65 6415 3188 **South Africa** 00 800 00 24 67 23 **Spain** 00 800 00 24 67 23 **Sweden** 00 800 00 24 67 23  
**Switzerland** 00 800 00 24 67 23 **Taiwan** 886 2 2578 7189 **Thailand** 66 2 651 8311 **United Arab Emirates** 36 1 459 6150 **United Kingdom** 00 800 00 24 67 23